



Product Code	28123
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2506591
CND	Z12159002
EAN/UPC	8023279281231
PARAF	989331804

## Description

### SMART NEBULIZER

Easy-to-use single button operated compressor nebulizer for home use. Light and compact design. Reliable pump with service life exceeding 1000 hours. It ensures highly efficient drug delivery for the treatment of asthma, allergies, and various respiratory diseases.

Supplied with inhalation tube, adult and child masks, inhalation mouthpiece, 2 air filters and multilingual manual in GB, FR, IT, ES, PT, DE, GR, RO, DK, HU, HR, BG.

## Technical Specifications

Max. capacity: 10 ml  
 Air flow: > 7 l/min  
 Operating mode: 20 min ON, 40 min OFF  
 Max. pressure: 210~400 Kpa  
 Nebulizing pressure, rate 60~180 Kpa; approx. 0.35 ml/min  
 Particle size: = 5µm 60%  
 Altitude operating range: 0-2000 m Noise level: approx. 55 dB  
 Voltage: AC 230 V, 50 Hz, 0.7A; = 80W  
 Size - weight: 140x140x100 mm - 1.1 kg  
 Safety: Class II, Type BF