



Product Code	28125
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2560193
CND	Z120105
EAN/UPC	6946594909333

Description

GIMA MANUAL SUCTION PUMP with adult/child catheter

Lightweight and compact aspirator to remove fluids from respiratory tracts or from the respiration support system and infective material from wounds.

- one hand operation
- reusable pump, disposable kit
- built-in integrated protection
- adjustable vacuum level for adult and child
- HEPA type microbial filter prevents contamination of the pump.

Class I Medical Device according to Regulation 2017/745/UE.

Supplied with disposable canister, with filter child and adult suction tube, adaptor and multilingual manual (GB, FR, IT, ES, PT, DE, PL, GR, RO, HU, HR).

Technical Specifications

Max vacuum: 450 mmHg
 Material: PE , PVC
 Airflow: >20 l/min
 Size: 180x220x80 mm
 Container capacity: 300 ml
 Weight: 275 g