



Product Code	28132
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1757691
CND	Z12159002
EAN/UPC	8033717011443
PARAF	926977428
GMDN	12716

## Description

HAPPYNEB II NEBULIZER - piston 230V- Faster-Jet

Happyneb stylish and portable piston nebulizer for a quick and efficacious aerosol therapy. Reliable, practical and safe. Noiseless and lubrication free. Unit includes an air filter to be replaced every 70 applications, and handle for easy carrying. It is built in accordance with European Standard EN 60601-1.

5 languages box and manuals: GB, FR, IT, ES, DE

Made in Italy.

## Technical Specifications

Capacity: 12 ml  
 Respirable volume < 5 µm 75%  
 MMAD (Mass Media Aerodynamic Diameter) 2.32 µm  
 Neb rate 0.45 ml/min  
 Residual volume 0.7 ml

## Standard accessories

Adult and pediatric face mask  
 Mouthpiece  
 Nasal prong  
 Fasterjet bulb  
 Connection tubing