

HOSPYNEB PROFESSIONAL NEBULIZER - piston - 230V -50/60 Hz



Product Code	28134
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1757691
CND	Z12159002
EAN/UPC	8033717011429
PARAF	903111868
GMDN	12716

Description

Hospyneb professional

Piston compressor nebulizers for professional and home care use. They embody all the dictates of modern aerosol therapy techniques. Reliable and simple to use, also by inexperienced users. Manual: GB, FR, IT, DE, ES.

Made in Italy.

Technical Specifications

- Operating voltage: 230 V - 50 Hz
- Size: 28 x 19 x h 17 cm
- Max pressure: 350 kPa (3.5 bar)
- Noise level: 60 dBA
- Operating pressure: 0-130 kPa (0-1.3 bar)
- Weight: 4.1 kg
- Air Flow: 15 l/min
- MDD class: IIa

Standard accessories

Fasterjet nebulizer with valve system
Nosepiece irrigation (polypropylene)
Mouthpiece (polypropylene) with valve (silicone)
Adult and child mask (PVC DEPH free)
Air tube (PVC DEPH free) and spare filter.