



Product Code	28135
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1757691
CND	Z12159002
EAN/UPC	8033717011450
PARAF	926977430
GMDN	12716

Description

HAPPYNEB III NEBULIZER - piston 230V

Happyneb stylish and portable piston nebulizer for a quick and efficacious aerosol therapy. Reliable, practical and safe. Noiseless and lubrication free. Unit includes an air filter to be replaced every 70 applications, and handle for easy carrying. It is built in accordance with European Standard EN 60601-1.

5 languages box and manuals: GB, FR, IT, ES, DE

Made in Italy.

Technical Specifications

Voltage: 230V/50 Hz
 Maximum pressure: 250 kPa
 Operating pressure: 90 kPa
 Air flow: 14 l/min
 Dimensions: 31x22x13 cm
 Weight: 2 kg
 Noise level 58 db
 Temporary use: Continuous use 30/30 min ON/OFF
 Class of risk - norms Ila - Class II, Type B

NEBULIZER TECHNICAL DATA

Capacity: 12 ml
 Respirable volume < 5 µm 75%
 MMAD (Mass Media Aerodynamic Diameter) 2.32 µm
 Neb rate 0.45 ml/min
 Residual volume 0.7 ml

Standard accessories

Nebjet bulb
 Adult face mask
 Pediatric face mask
 Nasal prong
 Mouthpiece
 Connection tubing