



Product Code	28141
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1379469
CND	A0680
EAN/UPC	8004007005367
PARAF	987474804
GMDN	58116

Description

GELID COAGULATIVE POWDER 50 g - box of 25

Super absorbent polymer powder used to turn infected medical liquid waste into a gel-like substance.

Mixed with liquid waste, it solidifies and reduces the risk associated with spills of organic liquids.

50 g of product have an absorption capacity of 3000 ml.

It can be used for surgical and bronchial suction system (canister or bags); canister for laboratory's liquid; dialysis liquid...

Material: sodium polycrylate cross-linked.

Packed in aluminum pouch.

Multilingual instructions: GB, FR, IT, ES, PT, FI, SK, PL, LV, RO, GR