



Product Code	28148
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	346438
CND	Z12080303
EAN/UPC	8033622540526
PARAF	903922843

Description

MAMILAT - BREAST PUMP - 230 V - 50 Hz

Easily transportable from one maternity ward to another, or for use at home. Made of highly heat-resistant, electrically insulated plastic material in conformity with the latest European safety standards. Working time: 20 minutes ON, 40 OFF. Suction adjustment knob. Supplied with 2 polycarbonate milk bottles (150 cc and 250 cc) antibacterial and hydrofobic filter, set of tubes 6x10 mm, breast shield, PVC air tube.

5 languages manual: IT, GB, ES, FR, DE.

Made in Italy.

Technical Specifications

Regulator for suction intensity: max suction -0.33 bar (-248 mmHg)
 Operating voltage: 220-230 V - 50/60 Hz - available on request 110 V - 60 Hz
 Max flow (air litres/min): 14 l/min
 Power consumption: 184 VA Noise level: 63 dB
 Size: 235x190x165 mm Weight: 2.2 kg
 Norms: CEI 62-5 (IEC601-1). IP Protection Class: IP21