



Product Code	28159
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1294815
CND	Z120105
EAN/UPC	8054619101414

Description

CLINIC PLUS Basic MPR

- Voltage: 230 V
- Autoclavable Jars: 2x5 l
- Footswitch: NO
- Flow direction diverter (Allows to direct suctioned liquids to any of the 2 jars): NO

Designed for professional aspiration of bodily fluids, tissues or bones of patients during or after surgery. Oilless and maintenance free piston type pumps provide high performances. Excellent suction capacities and max vacuum built up within a few seconds.

Available with vacuum gauge and autoclavable jars made of Makrolon® with 500 ml graduation.

High standards of safety in overflow protection system
Double security valve integrated in the jar and hydrophobic filter, safety trap bottle.

Main applications: EMERGENCY DEPT. / GENERAL SURGERY / GYNAECOLOGY / NEUROSURGERY / OBSTETRICS / OPERATING THEATRE / DENTAL PRACTICE

MPR (Multi Purpose Rail) models

This system enhances the versatility for easy and quick exchange of different accessories, with no need for tools. Equipped with five connections for rings of various diameters to fit jars of different sizes and types (2l, 3l, 5l), cannula holders and safety trap bottle.

Made in Italy.

Technical Specifications

Operating voltage: 230 V - 50/60 Hz
Maximum suction: - 0,90 bar (675 mmHg)
Operating cycle: continuous
Flow: 60 l/min
Power: 230 VA
Noise level: 51.7 dBA
Norms: CEI 62-5 (IEC 601-1); 93/42 EEC

Standard accessories

121°C autoclavable jar with overflow valve system: 2x5 l
Rings to accommodate 5 l jars: 2
Safety Trap Bottle (220 ml): 1
Antibacterial & Hydrophobic Filter (single-patient): 1
Autoclavable silicone tubes ø 8x14 mm length 150 cm: 1
Conical Connector ø 10-11-12mm: 1

Air suction inlet: 1

Footswitch with intermittent or continuous operation: NO

Change-Over System from jar to jar by soft-touch keys: NO

Power Cord with Schuko plug: 1

Multilingual manual: GB, FR, IT, ES ,DE, GR, PL, RO