

TRULON NON ABSORB. SUTURE gauge 5/0 circle 3/8
needle 19 mm - 45 cm - black



Product Code	22816
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	1781308
CND	H010201010102
EAN/UPC	8023279228168
GMDN	16584

Description

TRULON - polyamide nylon non absorbable monofilament suture
Equivalent Ethicon brand: ETHILON / NUROLON

- Gauge: 5/0
- Metric: 1.0
- Colour: Black
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 19 mm
- Suture length: 45 cm

Technical details:

- monofilament polyamide
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- for cuticular and subcuticular closure
- orthopedic surgery
- ophthalmic
- general closure
- micro surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.