

# YANKAUER CANNULA with bulb tip and suction tube 25 cm - sterile



Product Code	28174
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	46257
CND	A06010103
EAN/UPC	8004007005350
GMDN	47342

## Description

YANKAUER CANNULA with bulb tip and suction tube - 25 cm - box of 50

Ready-to-use rigid cannula for suction of fluids in surgical field.  
Anatomical and anti-slip handle for easy grip.  
Transparent to control the flow of suction fluids.  
Open tip with one central and four lateral holes (Ø 20 mm).  
Flexible, transparent and non-conductive suction tube.

Single use.  
Phthalates and latex free.

Packed in single pouch.  
Class IIa, sterile.