

# SPENCER JET COMPACT 300D PORTABLE SUCTION UNIT



Product Code	28177
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

RDM (NSIS)	79643
EAN/UPC	8023279281774

## Description

### JET COMPACT PORTABLE SUCTION UNIT

Extremely compact suction device to be used to remove fluids from the airway or respiratory support system.

The device creates a negative pressure that draws fluids through disposable tubing that is connected to the collection container.

Easy view of the setted vacuum thanks to the LED indicators.

Classification ISO: high vacuum / high flow.

Supplied with 12 V power cable, battery, disposable container with filter, disposable patient tube, manual GB, IT.

Made in Italy.

## Technical Specifications

Dimensions: 245 x 70 x 185 mm

IP grade: IP 12

Weight: 1.60 kg

Battery: 12 V DC Ni-MH

Vacuum range: from 50 to 550 mmHg

Battery operative time: 30 min ON, 30 min OFF

Flow (free flow): 32 l/min

Noise level: 70 dB

Battery recharging time: 5 h

Electrical connector requirement: 100-240 VAC 47-63 Hz 0.75 A max: 12 VDC

Container capacity: 300 ml

Patient tube length: 1 m