



Product Code	28188
Unit of sale	1 pc
Minimum order	1
Type	Spare part
Class	

## Description

Cover for jar 5 l  
- Autoclavable up to 143°C  
- Compatible with 28187

Autoclavable jars with overflow valve system

Equipped with screw cover with handle for easy grip, autoclavable silicone o-ring for tight seal, overflow valve system integrated in the cover, clear graduated scale in ml with 100 or 200 ml intervals.

Suitable for central vacuum systems and Gima suction units.

Available in two types of polycarbonate, Makrolon® and Apec®.