

# "SUPER VEGA" 36 BATTERY SUCTION UNIT



Product Code	28190
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1968723
CND	Z120105
EAN/UPC	8023279281903
PARAF	987472851
GTIN	08023279281903
GMDN	36777

## Description

SUPER VEGA BATTERY - plastic case - 36 l/min.

Professional, portable suction units for body, oral, nasal or tracheal liquids aspiration in adults or children. Supplied with three different operation options: AC/DC adapter, internal lead rechargeable battery and 12 V make it a versatile unit approved for NON-STOP operation without overheating. Visual and acoustic alarm indicating low battery level. Ergonomic handle ensures a comfortable and secure grip. Bottle with safety float control valve to prevent overflow.

Made in Italy.

## Technical Specifications

Suction power: 36 air l/min  
Max suction: -0.80 bar (600 mm/Hg)  
Operating Voltage: 100-240 V - 50/60 Hz  
Battery autonomy: approx 60 min  
Recharging time: 240 minutes  
Noise level: 65.5 dBA  
Size - Weight: 350x210x180 mm - 4.05 kg  
Norms: Class II Type BF. EN 60601-1-1. EN 60601-1-2. EN 60601-1-11. 93/42 EEC

## Standard accessories

Bottle 1 l - autoclavable at 121°C  
AC/DC universal adapter (100-240V - 50-60Hz)  
Silicone tubing set  
Anti-bacterial filter  
Suction catheter  
Cigarette lighter adaptor (12V 4A)  
User manual (GB, IT, FR, DE, ES, GR, PL, RO, PT)