

"SUPER VEGA" 36 BATTERY SUCTION UNIT



Product Code	28190
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1968723
CND	Z120105
EAN/UPC	8023279281903
PARAF	987472851
GMDN	36777

Description

SUPER VEGA BATTERY - plastic case - 36 l/min.

Professional, portable suction units for body, oral, nasal or tracheal liquids aspiration in adults or children. Supplied with three different operation options: AC/DC adapter, internal lead rechargeable battery and 12 V make it a versatile unit approved for NON-STOP operation without overheating. Visual and acoustic alarm indicating low battery level. Ergonomic handle ensures a comfortable and secure grip. Bottle with safety float control valve to prevent overflow.

Made in Italy.

Technical Specifications

Suction power: 36 air l/min
Max suction: -0.80 bar (600 mm/Hg)
Operating Voltage: 100-240 V - 50/60 Hz
Battery autonomy: approx 60 min
Recharging time: 240 minutes
Noise level: 65.5 dBA
Size - Weight: 350x210x180 mm - 4.05 kg
Norms: Class II Type BF, EN 60601-1-1, EN 60601-1-2, EN 60601-1-11, 93/42 EEC

Standard accessories

Bottle 1 l - autoclavable at 121°C
AC/DC universal adapter (100-240V - 50-60Hz)
Silicone tubing set
Anti-bacterial filter
Suction catheter
Cigarette lighter adaptor (12V 4A)
User manual (GB, IT, FR, DE, ES, GR, PL, RO, PT)