



Product Code	28206
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717011276
PARAF	981518564
GMDN	36777

Description

MINIASPEED BATTERY PRO - metal case - 30 l/min

High vacuum portable medical device for professional surgical aspiration.

- it can operate at 12 V (car) power supply, rechargeable batteries or directly by mains supply
- its two cylinders electro-compressor ensure high performances and long life
- up to 45 min. of autonomy with fully charged battery
- with 99% antibacterial hydrophobic filter
- equipped with autoclavable jar with safety valve, disposable suction liner and carrying bag.

Made in Italy.

Technical Specifications

Operating Voltage: 100-240 V 50/60 Hz AC - 14 V DC 4.28 A

Size: 370x100x240 mm

Weight: 5.8 kg (include kit)

Max suction: -0.85 bar

Suction power: 30 air l/min

Noise level: 65 dBA

Norms: MDD 93/42/ECC (class IIa - annex IX) IEC 60601-1, UNI ENISO 10079-1

Standard accessories

- Autoclavable jar 1000 cc with safety valve
- Silicone tubing set Ø 6x12 mm
- Disposable suction liner
- Sterile cannula
- Anti-bacterial hydrophobic filter
- Multivoltage battery charger
- Cigarette lighter power cord.
- Shoulder carrying bag
- User manual (GB and IT)