

# VEGA UNO SUCTION UNIT 1 liter jar GB,FR,IT,ES,PT,DE,CZ,LV,EE,FI,SI,BG



Product Code	28209
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	590729
CND	Z120105
EAN/UPC	8023279282092
PARAF	913758734
GTIN	08023279282092
GMDN	36777

## Description

### VEGA SUCTION UNIT

Model: Vega Uno  
Operating Power: 220-230 V - 50/60 Hz  
Flow: 15 l/m  
Manual: GB, FR, IT, ES, PT, DE, CZ, LV, EE, FI, SI, BG

A line of stylish, compact, and reliable aspirators for small surgery. Silent, high-flow, and maintenance-free pumps allow long product life and durability. Recommended for nasal aspiration, oral aspiration, and tracheal aspiration of body fluids (mucus, catarrh, or blood) in adults or children.

Equipped with continuously adjustable vacuum regulator and vacuum indicator  
Unbreakable 1,000 ml Makrolon® jar, autoclavable at 121°C  
Safety float control valve to prevent overflow  
Optional APEC® jar autoclavable at 143°C  
Manufactured with an electric and thermally insulated plastic body

Flovac® autoclavable jars and disposable collection system: in addition to the standard autoclavable jar, a disposable collection system is available, composed of a rigid reusable jar, a polyethylene disposable liner with hermetically sealed cover, and a hydrophobic, antireflux, antibacterial filter operating as an overflow valve.

Made in Italy.

## Technical Specifications

Operating Voltage: 220-230 V - 50/60 Hz, 120-127 V - 60 Hz (on request)  
Power Consumption: 184 VA  
Bottle Capacity: 1,000 ml  
Flow (air litres/min): 15 l/min  
Max Suction: -0.75 bar (563 mm/Hg)  
Operating Cycle: Continuous  
Noise Level: 59.6 dBA  
Size - Weight: 25x19x16.5 cm - 2.2 kg  
Norms: EN 60601-1-1, EN 60601-1-2, EN 60601-1-11, 93/42/EEC

## Standard accessories

Autoclavable jar (1 l) with overflow valve system  
Antibacterial and hydrophobic filter  
Set of silicone tubes 6x10 mm  
Conical connector for probes