

TRULON NON ABSORB. SUTURE gauge 6/0 circle 3/8  
needle 16 mm - 45 cm - black



|               |                |
|---------------|----------------|
| Product Code  | 22822          |
| Unit of sale  | box of 12 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | III            |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 1781308       |
| CND        | H010201010102 |
| EAN/UPC    | 8023279228229 |
| GMDN       | 16584         |

Description

TRULON - polyamide nylon non absorbable monofilament suture  
Equivalent Ethicon brand: ETHILON / NUROLON

- Gauge: 6/0
- Metric: 0.7
- Colour: Black
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 16 mm
- Suture length: 45 cm

Technical details:

- monofilament polyamide
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- for cuticular and subcuticular closure
- orthopedic surgery
- ophthalmic
- general closure
- micro surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.