

# COVER+VALVE for bottle 5 litres (28221)



Product Code	28225
Unit of sale	1 pc
Minimum order	1
Type	Spare part
Class	

EAN/UPC 8023279282252

## Description

Cover for jar 5 l  
- Autoclavable up to 121°C  
- Compatible with 28221

Autoclavable jars with overflow valve system  
Equipped with screw cover with handle for easy grip, autoclavable silicone o-ring for tight seal, overflow valve system integrated in the cover, clear graduated scale in ml with 100 or 200 ml intervals.  
Suitable for central vacuum systems and Gima suction units.  
Available in two types of polycarbonate, Makrolon® and Apec®.