



Product Code	28245
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717010767
PARAF	903150098
GMDN	36777

Description

Piston-type continuous cycle electric aspirators give high performance and great durability. Equipped with a protective thermal cut-out relay. They require no maintenance or lubrication. A motor-protector cap totally prevents aspirated liquids or secretions from reaching and damaging the vacuum pump. Ideal for clinical use. Made in Italy.

- Aspeed 22 l
- Power: 230 V
- Pump: double
- Case: metal

Technical Specifications

Operating voltage: 230 V-50 Hz
 Bottle capacity: 1 l
 High vacuum: low flow
 Air flow: 22 l/min
 Adjustable vacuum level: 0÷ -0.85 bar (0÷ -85 kPa)
 Weight: 4.2 kg
 Case material: metal
 Noise level: 63 dBA

Standard accessories

Autoclavable polycarbonate jar 1,000 cc with safety valve (overflow protection)
 Disposable suction liner 1 l
 Antibacterial hydrophobic filter
 Sterile disposable cannula
 Sterile manual flow regulator
 Set of atoxic sterilizable silicone tubes
 Power Cable
 User Manual: GB, FR, IT, DE, ES