

ADAPTOR FOR CATHETER for ASPEED suction aspirator



| | |
|---------------|----------------|
| Product Code | 28257 |
| Unit of sale | box of 100 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I S |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 974537 |
| CND | A0701 |
| EAN/UPC | 4034342046040 |

Description

Adaptor to connect catheter to silicone tube for ASPEED

- Box of 100 pcs.