

# "ASPEED 3" SUCTION ASPIRATOR - 230V single pump - 1 l



Product Code	28279
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717011405

## Description

ASPEED 3.0 aspirator 16 l/min - 1 l

Portable, compact reliable aspirator ideal for small surgery.  
Piston type continuous cycle electric medical device.  
High vacuum low flow.  
No maintenance or lubrication.  
Made in Italy.

Can be used also with optional 2-liter jar (code 28294).

## Technical Specifications

Operating voltage: 230 V - 50 Hz  
Power consumption: max 150 VA  
Bottle capacity: 1 l  
High vacuum: low flow  
Air flow: 16 l/min  
Adjustable vacuum level: 0÷ -0.65 bar (0÷ -65 kPa)  
Size: 300x170xh 260 mm  
Weight: 2.3 kg  
Noise level: 65 dBA  
Operating cycle: continuous  
Class: II  
Type: BF

## Standard accessories

Autoclavable polycarbonate jar 1,000 cc or 2,000 cc with safety valve (28258 or 28294)  
Disposable liner 1 l (28259) only for 28258  
Antibacterial hydrophobic filter  
Sterile disposable cannula  
Sterile manual flow regulator  
Set of atoxic sterilizable silicone tubes  
Power Cable  
User Manual: GB. FR. IT. DE. ES