



Product Code	28280
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717011382
PARAF	987471240

## Description

Piston-type continuous cycle electric aspirators give high performance and great durability. Equipped with a protective thermal cut-out relay. They require no maintenance or lubrication. A motor-protector cap totally prevents aspirated liquids or secretions from reaching and damaging the vacuum pump. Ideal for clinical use. Made in Italy.

- Aspeed 2 15 l
- Power: 230 V
- Pump: single
- Case: plastic

## Technical Specifications

Operating voltage: 230 V-50 Hz  
 Bottle capacity: 1 l  
 High vacuum: low flow  
 Air flow: 15 l/min  
 Adjustable vacuum level: 0÷ -0.85 bar (0÷ -85 kPa)  
 Weight: 2,7 kg  
 Case material: plastic  
 Noise level: 63 dBA

## Standard accessories

Autoclavable polycarbonate jar 1,000 cc with safety valve (overflow protection)  
 Disposable suction liner 1 l  
 Antibacterial hydrophobic filter  
 Sterile disposable cannula  
 Sterile manual flow regulator  
 Set of atoxic sterilizable silicone tubes  
 Power Cable  
 User Manual: GB, FR, IT, DE, ES