



Product Code	28282
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717011283

## Description

### MINIASPEED EVO BATTERY SUCTION - 1 I

Professional portable high flow - high vacuum aspirators with integrated rechargeable batteries.  
It works at 12 V power supply with batteries or directly by mains supply, continuously at 110-230 V while recharging batteries.  
Up to 45 minutes of autonomy with fully charged battery.

Made in Italy.

## Technical Specifications

Operating voltage: with internal battery DC 12 V - 4 Ah with power supply 100÷240 V ~ 50/60 Hz - 14 VDC 4,28A

Power consumption: max 50 W

Battery autonomy: about 45 minutes, recharging time approx 6 hours

Size: 134x378xh 263 mm

Weight: 4.7 kg

Bottle capacity: 1,000 ml

Flow: 30 air litres/min

Noise level: 65 dBA at 1 m

Operating cycles: continuous (45 min battery autonomy)

MAX vacuum: 0 ± 85 kPa

Norms: MDD 93/42/EEC (Class IIa - annex IX) EN60601-1(IEC 60601-1) / EN60601-1-2 (IEC 60601-1-2)

UNI EN ISO 10079-1/ UNI CEI EN ISO 14971

Degree of protections against ingress of liquids: IP21

Class IIa - Type BF

## Standard accessories

Autoclavable 1 l jar with safety valve

2 bacterial hydrophobic filters

Sterile disposable cannula CH14

User manual: GB, IT, FR, DE, ES

Sterile disposable manual regulator

Power supply cord - 2 m

Disposable silicone liner 1 l

Shoulder carrying bag

Multivoltage battery charger

Cigarette lighter power cord

Silicone tubing set Ø 6x12 mm