

MINI ASPEED EVO BATTERY SUCTION 1 I with ambulance bracket



Product Code	28283
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	1758469
CND	Z120105
EAN/UPC	8033717011306

Description

MINIASPEED EVO BATTERY SUCTION - 1 I for ambulance

Professional portable high flow - high vacuum aspirators with integrated rechargeable batteries. It works at 12 V power supply with batteries or directly by mains supply, continuously at 110-230 V while recharging batteries. Up to 45 minutes of autonomy with fully charged battery.

Made in Italy.

Technical Specifications

Operating voltage: with internal battery DC 12 V - 4 Ah with power supply 100÷240 V ~ 50/60 Hz - 14 VDC 4,28A
Power consumption: max 50 W
Battery autonomy: about 45 minutes, recharging time approx 6 hours
Size: 134x378xh 263 mm
Weight: 5 kg
Bottle capacity: 1,000 ml
Flow: 30 air litres/min
Noise level: 65 dBA at 1 m
Operating cycles: continuous (45 min battery autonomy)
MAX vacuum: 0 ± 85 kPa
Norms: MDD 93/42/EEC (Class IIa - annex IX) EN60601-1(IEC 60601-1) / EN60601-1-2 (IEC 60601-1-2)
UNI EN ISO 10079-1/ UNI CEI EN ISO 14971
Degree of protections against ingress of liquids: IP21
Class IIa - Type BF

Standard accessories

Autoclavable 1 l jar with safety valve
2 bacterial hydrophobic filters
Sterile disposable cannula CH14
User manual: GB, IT, FR, DE, ES
Sterile disposable manual regulator
Power supply cord - 2 m
Disposable silicone liner 1 l
Shoulder carrying bag
Wall bracket for ambulance
Multivoltage battery charger
Silicone tubing set Ø 6x12 mm