

# MAXI ASPEED SUCTION 60 I 2x4 I jar with footswitch - 230V



Product Code	28288
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717011344

## Description

### MAXI ASPEED SURGICAL ASPIRATORS

- Maxi Aspeed 60 I
- Jar: 2x4 l
- Power: 230 V
- Foot switch: YES
- Flow direction diverter (allows to direct suctioned liquids to any of the 2 jars): YES

Surgical aspirators for aspiration of body liquids in operating theatres and hospitals.

Easy to carry  
Equipped with 4 antistatic wheels with brake and carrying telescopic handle.

Developed to run continuously  
Made in highly heat resistant, electrically insulated plastic material in conformity with latest European safety standards.

Vacuum regulator with BAR/PSI vacuumeter and automatic flow direction regulator to easily switch from a jar to the other.

Supplied with 2 autoclavable polycarbonate jars with safety valve (overflow protection), set of atoxic sterilizable silicone tubes, antibacterial hydrophobic filters (2) and multilingual manual (GB, FR, IT, DE, ES).

Made in Italy.

## Technical Specifications

- Operating voltage: 230 V - 50/60 Hz (Available on request 110 V-60 Hz)
- Air flow: 60 L/min
- Operating cycle: continuous
- Adjustable vacuum level: 0-0.90 bar (0-90 kPa)
- Size: 470 x 580 x h 560 mm
- Weight: 15 kg

## Standard accessories

- Foot switch
- Flow direction regulator