



Product Code	28289
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717011351

## Description

### MAXI ASPEED SURGICAL ASPIRATORS

- Maxi Aspeed 90 I
- Jar: 2x2 l
- Power: 230 V
- Foot switch: NO
- Flow direction diverter (allows to direct suctioned liquids to any of the 2 jars): YES

Surgical aspirators for aspiration of body liquids in operating theatres and hospitals.

Easy to carry  
Equipped with 4 antistatic wheels with brake and carrying telescopic handle.

Developed to run continuously  
Made in highly heat resistant, electrically insulated plastic material in conformity with latest European safety standards.

Vacuum regulator with BAR/PSI vacuometer and automatic flow direction regulator to easily switch from a jar to the other.

Supplied with 2 autoclavable polycarbonate jars with safety valve (overflow protection), set of atoxic sterilizable silicone tubes, antibacterial hydrophobic filters (2) and multilingual manual (GB, FR, IT, DE, ES).

Made in Italy.

## Technical Specifications

Operating voltage: 230 V - 50/60 Hz  
 Noise level: 50-60 dBA  
 Air flow: 90 l/min  
 Operating cycle: continuous  
 Adjustable vacuum level: 0-0.90 bar (0-90 kPa)  
 Size: 470x580x560 mm  
 Weight: 17 kg