



Product Code	28296
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No
RDM (NSIS)	387280
CND	Z120606
EAN/UPC	8019781472205
UMDNS	12415

## Description

MAG 700 MAGNETOTHERAPY - 1 channel - 11 programmes

New portable magnetotherapy ideal for home therapy. MAG 700 offers a high level of professional performances and it's indicated for several pathologies. It allows comfortable home therapy through its basic accessories. MAG 700 has 7 pre-set programmes for the treatment of several pathologies, 3 programmes identified by the working frequency (from 5 to 20 Hertz) and an Autoscan programme with continuous frequency variation for the treatment of both hard and soft tissues at the same time. Possibility to adjust the magnetic field intensity for all programmes up to 70 Gauss. Working frequencies are between 5 and 100 Hertz.

Supplied with elastic therapeutic belt with 2 solenoids, power supply, carrying bag and manual (GB, FR, IT, ES). On request DE.

## Technical Specifications

11 PROGRAMMES

7 pre-set programmes:

- Osteoarthritis
- Osteoporosis
- Arthritis
- Fractures
- Articular pain
- Muscle pain
- Anti-inflammatory

3 working freq.:

- 5 Hz
- 10 Hz
- 20 H

AUTOSCAN