



Product Code	28308
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A
RDM (NSIS)	1995066
CND	II A
EAN/UPC	8019781413031

## Description

POWERSONIC BATTERY ULTRASOUND THERAPY

The Powerersonic ultrasound-therapy is a portable device ideal for treating conditions and diseases affecting the musculoskeletal system, providing analgesic and antiinflammatory relief across a wide range of applications. The device is suitable for use by both patients and healthcare professionals.

Supplied with medical power supply, rechargeable battery, ultrasound gel, carrying and storage bag, user manual in GB, IT, FR, ES, DE on request.

## Technical Specifications

Programmes: 13 programmes  
 Treatment time: 5, 10, 15 minutes  
 Ultrasound frequency: 1MHz ultrasound  
 Maximum output power: 6.4W (with power supply) and 3.2W (with battery)  
 Maximum intensity: 1.6W/cm<sup>2</sup> (with power supply) and 0.8W/cm<sup>2</sup> (with battery)  
 Beam type: collimated  
 3 ultrasound beam intensity levels: L (low), M (medium) e H (high)  
 Power supply: 100-240V, 50/60 Hz or Ni-MH AAA 850mAh 4.8 V battery  
 Size: 204x63x58 mm  
 Protection: class II type BF applied part, IP22, IPX7  
 Direct contact device (not suitable for water immersion)

**TREATMENTS**

- Myofascial pain
- Low back pain
- Tendinitis
- Muscle pain
- Arthrosis
- Cubital tunnel syndrome
- Carpal tunnel syndrome
- Epicondylitis (Tennis elbow)
- Lumbar stenosis
- Venous and pressure ulcers
- Epitrochleitis (Golfer's elbow)
- Sciatica