



Product Code	28309
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2586237
CND	Z120610
EAN/UPC	8019781472175
UMDNS	11248
PARAF	921453294

Description

MIO-SONIC ULTRASOUND THERAPY

Innovative and high performance device both for home use and professional use.

It combines compactness and high performances in a small, light and easy to use applicator. The lightness and handiness allows a continuous use without any difficulty.

Mio Sonic can be used for the following applications:

- pain therapy
- sport injuries
- aesthetic treatments

Supplied with gel bottle, carrying and storage bag, medical power supply and multilanguage manual: GB, FR, IT, ES. On request DE.

Technical Specifications

9 programs;

1 MHz ultrasound, maximum output power 6.4 W, maximum intensity 1.6 W/cm²;

Collimated ultrasound beam for an optimal penetration beam;

3 ultrasound beam intensity levels: 0,08 W/cm² (low) - 0,8 W/cm² (medium) - 1,6 W/cm² (high);

Power supply: 100-240V, 50/60 Hz

Direct contact device (not suitable for water immersion)

Size-weight: 200x47x69 mm - 190 g