


GELLED ELECTRODES 45x80 mm



Product Code	28317
Unit of sale	box of 4 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2269235
CND	N010299
EAN/UPC	8033695723581
PARAF	987472848
GMDN	17191

Description

DISPOSABLE GELLED ELECTRODES - latex free

High material flexibility allows easy application of electrodes on treated areas. Solid gel maintains the electrode adhesion and conductivity unchanged during procedure.