



Product Code	28318
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2306822
CND	Z12069099
EAN/UPC	8019781112507

Description

OSTEOMAT MATTRESS FOR MAGNETOTHERAPY

Osteomat is composed of 4 active sectors, containing 2 solenoids each, for a total of 8 professional solenoids. The Osteomat is the ideal solution for total body treatment in a single session, thanks to easy-to-use mode and extreme flexibility. The configuration of Osteomat stimulates the body areas (shoulders, column, pelvis) that are more affected by bone mass reduction. This is typical of osteoporosis pathology. Suitable for MAG 2000, MAG 2000 PLUS, MAG 3000 and LaMagneto models. Supplied with non-woven fabric cover (pack of 2 pieces).

Suggested for night therapy during sleep, day therapy in paralyzed patients, active patients that don't tolerate long treatments with therapeutic belt and prefer night therapy.

4 active sectors (2 solenoids for each sector, totally 8 solenoids)

Manual: GB, FR, IT, ES