



Product Code	28319
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717010507

Description

ASPEED 3.0 aspirator 16 l/min - 2 l

Portable, compact reliable aspirator ideal for small surgery.
Piston type continuous cycle electric medical device.
High vacuum low flow.
No maintenance or lubrication.
Made in Italy.

Technical Specifications

Operating voltage: 230 V - 50 Hz
Power consumption: max 150 VA
Bottle capacity: 2 l
High vacuum: low flow
Air flow: 16 l/min
Adjustable vacuum level: 0÷ -0.65 bar (0÷ -65 kPa)
Size: 300x170xh 260 mm
Weight: 2.3 kg
Noise level: 65 dBA
Operating cycle: continuous
Class: II
Type: BF

Standard accessories

Autoclavable polycarbonate jar 1,000 cc or 2,000 cc with safety valve (28258 or 28294)
Disposable liner 1 l (28259) only for 28258
Antibacterial hydrophobic filter
Sterile disposable cannula
Sterile manual flow regulator
Set of atoxic sterilizable silicone tubes
Power Cable
User Manual: GB, FR, IT, DE, ES