

TRUSILK NON ABSORB. SUTURE gauge 2/0 circle 3/8 needle 26 mm - 45 cm - black



Product Code	22832
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No
RDM (NSIS)	1781309
CND	H010202020302
EAN/UPC	8023279228328
PARAF	952038863
GMDN	13910

Description

TRUSILK - braided black silk non absorbable multifilament suture
Equivalent Ethicon brand: MERSILK / PERMA-HAND SILK

- Gauge: 2/0
- Metric: 3.0
- Colour: Black
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 26 mm
- Suture length: 45 cm

Technical details:

- made from 100% protein fiber spun
- coating of wax/silicone
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- general surgery
- ophthalmic surgery
- cardiovascular surgery
- gastrointestinal surgery
- orthopedic surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.