

ELASTIC BELT 100x8 cm for reusable electrodes



Product Code	28323
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2432642
CND	C0280
EAN/UPC	8033003265024
PARAF	987472053
GMDN	35709

Description

UNIVERSAL ELASTIC BELTS FOR REUSABLE ELECTRODES

Elastic belts - blue - 100x8 cm

Elastic belts for reusable electrodes are easy to use and assure a very good seal.
Made of synthetic fibres (nylon, elastomer and transparent nylon yarn) to avoid skin irritation in patients.
Belts are available in 4 different lengths to match almost all kinds of reusable electrodes on the market.