

# ELASTIC BELT 60x8 cm for reusable electrodes



Product Code	28331
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2432640
CND	C0280
EAN/UPC	8033003264973
GTIN	08033003264973

## Description

### UNIVERSAL ELASTIC BELTS FOR REUSABLE ELECTRODES

Elastic belts - blue - 60x8 cm

Elastic belts for reusable electrodes are easy to use and assure a very good seal.  
Made of synthetic fibres (nylon, elastomer and transparent nylon yarn) to avoid skin irritation in patients.  
Belts are available in 4 different lengths to match almost all kinds of reusable electrodes on the market.