



Product Code	28338
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2586120
CND	Z120606
EAN/UPC	8019781472090
PARAF	921487132
GMDN	47593

## Description

MAG 1000 MAGNETOTHERAPY - 2 channels - 20 programmes

Mag 1000 is a high frequency PEMF device particularly recommended for the treatment of pathologies affecting soft tissues like muscles and tendons. It offers 20 pre-set programmes with 3 adjustable intensity levels. The device generates an electric signal with a carrier frequency of 27 MHz, which is modulated with specific frequencies (from 8 Hz to 640 Hz) and therefore is able to soothe the anatomical site affected by pain. It is projected and indicated for healing, rehabilitation and functional recovery treatments of many pathologies. Continuous operating mode. It is intended for both professional and homecare use.

Size: 153x91x38 mm

Supplied with elastic therapeutic belt, carrying bag, rechargeable battery pack 2,300mAh, battery charger, user manual in GB, FR, IT, ES, DE.

## Technical Specifications

20 programmes:  
 Arthritis-Arthrosis  
 Bruises  
 Rheumatisms  
 Myalgia  
 Pain-articular prosthesis  
 Tendonitis  
 Herniated disc  
 Hematoma  
 Muscular pain  
 Epicondylitis  
 Osteoporosis  
 Epytrocleitis  
 Fractures  
 Lumbago  
 Sprains  
 Muscular contractures  
 Dislocations 2h  
 Medium application  
 Traumas 4h  
 Long application