



Product Code	28338
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2586120
CND	Z120606
EAN/UPC	8019781472090
PARAF	990983140
GMDN	47593

Description

MAG 1000 MAGNETOTHERAPY - 2 channels - 20 programmes

Mag 1000 is a high frequency PEMF device particularly recommended for the treatment of pathologies affecting soft tissues like muscles and tendons. It offers 20 pre-set programmes with 3 adjustable intensity levels. The device generates an electric signal with a carrier frequency of 27 MHz, which is modulated with specific frequencies (from 8 Hz to 640 Hz) and therefore is able to soothe the anatomical site affected by pain. It is projected and indicated for healing, rehabilitation and functional recovery treatments of many pathologies. Continuous operating mode. It is intended for both professional and homecare use.

Size: 153x91x38 mm

Supplied with elastic therapeutic belt, carrying bag, rechargeable battery pack 2,300mAh, battery charger, user manual in GB, FR, IT, ES, DE.

Technical Specifications

20 programmes:
 Arthritis-Arthrosis
 Bruises
 Rheumatisms
 Myalgia
 Pain-articular prosthesis
 Tendonitis
 Herniated disc
 Hematoma
 Muscular pain
 Epicondylitis
 Osteoporosis
 Epytrocleitis
 Fractures
 Lumbago
 Sprains
 Muscular contractures
 Dislocations 2h
 Medium application
 Traumas 4h
 Long application