



Product Code	28342
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1225279
CND	Z120601
EAN/UPC	8023279283426
UMDNS	13-782
PARAF	977812270

Description

GIMA 9 TENS - 2 channels

Pocket TENS with 9 programmes and 2 channels. 9 adjustable intensities.
Relieves muscle pain and stimulates blood circulation.
Frequency: 2-80 Hz.

Supplied in a nylon bag with 2AAA batteries and 4 electrodes.

Size: 128 x 56 x 16 mm

LCD size: 48 x 25.4 mm

Weight: 67 g

Manual: GB, FR, IT, ES, DE, PT, GR, PL.