



|               |                   |
|---------------|-------------------|
| Product Code  | 28347             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279283471 |
| PARAF   | 950851079     |

## Description

### DIATER CONDUCTIVE CREAM 1 L - HYALURONIC ACID

Essential, simple and effective conductive creams. Rich in active ingredients which are recognized in musculoskeletal therapies, in support of the functional recovery of athletes and in the aesthetic treatments of epidermis and dermis. Exceptional absorption and therapies enhancement keep the activity focused on the portion of treated tissue, with marked stability and effectiveness.

Absence of parabens or petrolatum.

Neutral base conductive cream moisturizes and protects tissue; suitable for mixes and cocktails.

Arnica Montana conductive cream decongests and protects tissue.

For acute and chronic inflammation of tissues or joints and Lipodystrophy.

Hyaluronic acid conductive cream for hydro balancing of the tissues.

Protective of cartilage and against phlegm and tissue trauma spills.

Treatment of scars, hypotonia and post-surgery outcomes.

Suitable for passive mobilizations.

Creams are preserved into thermo/photo protected PET containers with cap and undercap and dispenser.

For professional use only in different fields: Physiotherapy, Physiatry, Plastic Surgeons, Aesthetics, Orthodermic therapy, Athletic Preparation, Massage and Podiatry.

Label in: GB, IT, FR, ES.

Made in Italy.