


# GELLED ELECTRODES 45x80 mm with cable



Product Code	28373
Unit of sale	box of 4 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2289174
CND	N010299
EAN/UPC	8033003260265
PARAF	984629131
GMDN	17191

## Description

DISPOSABLE GELLED ELECTRODES WITH CABLE - latex free

High material flexibility allows easy application of electrodes on treated areas. Solid gel maintains the electrode adhesion and conductivity unchanged during procedure.

The electrodes are equipped with 10 cm long cable, with 2 mm socket.