



Product Code	28382
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1225566
CND	Z120601
EAN/UPC	8019781501011
UMDNS	13762
PARAF	926559523

Description

I-TECH PHYSIO PROFESSIONAL STIMULATOR

Professional stimulator with 2 independent channels

- compensated biphasic square wave: current quantity fro positive pole to negative one is constant to avoid hazardous thermal effect of polarization
- monophasic square wave
- incontinence and ionophoresis programmes
- denervated treatment
- programmable TENS - EMS memories
- adjustable frequency, width impulse, time, ramp
- remaining therapy time and programmes instructions on the LC D
- maximum intensity 200 mApp;

Can use same gelled electrodes with cable as Mio-Care

Supplied in a carrying bag with:

- rechargeable battery and battery charger
 - 1 set of 4 pregelled electrodes 48x48 mm
 - 1 set of 4 pregelled electrodes 50x90 mm
 - 2 silicone conductive electrodes 60x90 mm
 - 2 sponges 70x100 mm
 - 1 elastic belt for electrodes fastening
 - 2 cables, 2 mm connection male plug
 - 4 cable splitter 2 mm male and female plug
- User manual: GB, FR, IT, ES, on request PT, DE.

Backlit display

Digital easy to use keyboard

Technical Specifications

PROGRAMMES:

- 21 NEMS
- 14 TENS
- 15 BEAUTY
- 9 UROLOGY
- 27 REHABILITATION
- 13 FREE MEMORIES