


GELLED ELECTRODES 45x35 mm with cable



Product Code	28384
Unit of sale	bag of 4 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2289167
CND	N010299
EAN/UPC	8033003260272
PARAF	977812282
GMDN	17191

Description

DISPOSABLE GELLED ELECTRODES WITH CABLE - latex free

High material flexibility allows easy application of electrodes on treated areas. Solid gel maintains the electrode adhesion and conductivity unchanged during procedure.

The electrodes are equipped with 10 cm long cable, with 2 mm socket.