

# CLIPS ADAPTORS 4 mm - red



|               |                  |
|---------------|------------------|
| Product Code  | 28393            |
| Unit of sale  | box of 10<br>pcs |
| Minimum order | 1                |
| Type          | Medical device   |
| Class         | I                |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2237241       |
| CND        | Z12050380     |
| EAN/UPC    | 8033003263020 |
| PARAF      | 977812344     |
| GMDN       | 40811         |

## Description

4 mm CLIP ADAPTORS - red - box of 10