



Product Code	28401
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1456070
CND	Z120601
EAN/UPC	8019781501080
UMDNS	13762
PARAF	936980580

Description

T-ONE COACH ELECTROTHERAPY - 4 channels

A high performant muscle stimulator designed to improve physical activity. Easy to use programs ideal device for not professional athletes and for sports lovers. It is not only a muscle stimulator equipment. It is supplied by TENS protocols for pain relief and injuries recovery.

Main features

- constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- four independent channels unit (up to 8 connectable electrodes)
- compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect;
- monophasic wave generator (only for iontophoresis programs)
- 2+2 users programs
- backlight graphic colour display
- user feedback: the device recognizes patient connection
- soft-touch easy to use keyboard
- on-line 5 languages guide (GB, FR, IT, ES, DE) on display
- 4,8 V rechargeable battery pack
- personalized programs for man and woman
- maximum intensity 240 mApp, with 1 mA step regulation

Size: 150x90x35 mm

Technical Specifications

• BEAUTY (18 PROGRAMMES)

- Firming up
- Toning up
- Definition
- Modelling
- Microlifting
- Lipolysis abdomen
- Lipolysis thighs
- Lipolysis glutei and hips
- Lipolysis arms
- Tissue elasticity
- Capillarization
- Heaviness in legs
- Face Microlifting Handy Trode

Breast Microlifting Handy Trode
Décolleté Microlifting Handy Trode
Face wrinkles Handy Trode
Forehead Wrinkles Handy Trode
Collagen

• TENS (16 PROGRAMMES)

Tens eapid
Tens Endorphinic
??Tens Max Values
Inflammatory
Neck pain
Backache/Sciatic pain
Sprains/Bruises
Vascularization
Muscle Relaxant
Hand and wrist pain
Plantar stimulation
Epicondylitis
Epitroclea
Periarthritis
Tens Burst
Microcurrent

• NEMS (14 PROGRAMMMES)

Warming up
Resistance
Resistance strength
Basic strength
Fast strength
Explosive strength
Deep capillarization
Muscle recovery
Agonist/Antagonist
Sequential tonic contractions
Sequential phasis contractions
Muscle Relaxant
Deep massage
EMS

• 2+2

Tens rapid
Tens endorphinic
firming up
Toning up
Definition
Modelling
Microlifting
Muscle Relaxant
Massage
Capillarization

Standard accessories

4 connection cables (8 connectable electrodes)
Kit of pre-gelled self-adhesive electrodes (square, rectangular, round)
Ionophoresis kit: elastic belt, 2 silicone electrodes, 2 sponges
Internal Ni-Mh rechargeable battery pack and battery charger
Electrode position manual and user manual (GB, FR, IT, ES). On request DE, PT