



Product Code	28406
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2230535
CND	Z120601
EAN/UPC	8023279284065
PARAF	986873469

## Description

GIMA-THER TENS - 36 medical programmes

Gima-Ther is a dual channel stimulator with 3 basic functions: TENS, EMS and MASSAGE.

It has 36 programmes (18 TENS, 15 EMS and 3 Massage) and applies electric currents in the low-frequency range for therapy. It helps with blood circulation and to relieve chronic and postoperative pains.

Recommended for muscle training and also for body building.

Electric impulses are transcutaneous and transmitted to nerves, muscle fibers or tissues through electrode. The intensity of the dual channel can be adjusted independently and can be applied individually to one body part.

Supplied with 4 pads 50x50 mm, 2 output cables, 4 AAA batteries and instruction manual (GB, FR, IT, ES, DE).

## Technical Specifications

### 36 PROGRAMMES - TENS/EMS/MASSAGE:

TENS (Transcutaneous Electrical Nerve Stimulation)

Body parts: shoulder, low back, joint, hip, arms, legs

Treating programme: 18 programmes

Treating time: 30 minutes

EMS (Electrical Muscle Stimulation)

Body parts: shoulder, back, limb and hip

Treating programme: 15 programmes

Treating time: 28, 30 and 32 minutes

MASSAGE

Treating programme: 3 programmes

Treating time: 30 minutes

Technical specifications:

Channel 2

Waveform: Bi-phase square-wave pulse

Output voltage: Max. 120 Vpp (at 1000 ohm load)

Output current: Max. 60 mA (at 1000 ohm load)

Output intensity: 0 to 60 levels, adjustable

Output characteristics: Constant Current (CC)

Power supply: 6.0 V D.C., 4x AAA batteries 1.5 V

Display: Large backlit LCD

Dimension - Weight: 120.5x69.5x27 mm - 104 g

