

# ANYCARE DIGITAL AIR COMPRESSION MASSAGE with 2 leg sleeve A



Product Code	28415
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8809315671514
PARAF	987476138

## Description

### ANYCARE DIGITAL COMPRESSION MASSAGE - LEG A

Pressure intensity and time can be adjusted according to user's needs and doctor's indications.

Individual chambers pressure testing.

Portable and easy operation. For hospital, clinics and homecare use.

Supplied with 2 legs A (74x62 cm) 4 chambers, 1 double hose and 2 plastic feet guides, multilingual manual (GB, FR, IT ,ES).

Power consumption: 32 W (40 VA) +10%

Voltage: AC 220V - 50-60 Hz

Human body has a peripheral circulation of blood and lymph. Air compression therapy is an automatic drainage and circulation system which helps patients to recover from their disease by improving blood and lymph circulation and metabolism. It consists of air chambers and an air pump offering air compression massage to legs, arms or waist.

Home care indications:

- Cold feet or hands
- Settlement of stress
- Exercise effects
- Insomnia
- Healthier body and elastic skin
- Old age, weakness and fatigue