

# KURA TENS - 2 channels - 6 programmes IT,DK,LT,FI,HU,RO,HR,SI,GR,BG



Product Code	28432
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2663092
CND	Z120601
EAN/UPC	8023279284324
PARAF	951163411

## Description

KURA TENS - 6 programmes - IT, DK, LT, FI, HU, RO, HR, SI, GR, BG

Dual channel electrotherapeutic device with massage feature ideal for pain relief associated with aching muscles due to physical exertion or regular work and domestic activities. It is also aimed at providing symptomatic relief and managing chronic and intractable pain, as well as pain relief related to arthritis.

- 2 channels with 4 electrode pads
- 15 intensities per channel
- 6 main modes (composite, thump, press, slap, acupuncture, muscular relaxation)
- 24 sub-modes (4 for each main mode)
- large LCD with blue backlight.

Supplied with 2 wires (120 cm), 4 electrode pads (50x50 mm) and multilingual manual (see above).

## Technical Specifications

Output channels: 2  
Body parts: waist, back, arm, leg, joint, shoulder, adductor and gluteus, chest, belly, feet  
Treatment duration: 15 to 30 minutes  
Band width: 0~100 us  
Pulse frequency range: 0-100 Hz  
Output voltage: 0~60 V ± 20% (500 ohm load)  
Power supply: 4 x 1.5 V AAA batteries  
Size - Weight: 120.3x60.3x20.6 mm - 73 g (batteries excluded)  
Classification: internally powered, Type BF applied part, IP22