

DOCTOR LIFE MK400 PROFESSIONAL COMPRESSION SYSTEM with 2 legs



Product Code	28441
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1379293
CND	Z120607
EAN/UPC	8809802641235
GMDN	10969

Description

MK400L professional compression system with 2 XL legs, and 2 single hoses.

MK400L is a professional air compression therapy system with essential functions for different kinds of patients. This system provides 3 compression modes (wave, squeezing, double squeezing), reverse compression settings, skip functions to select directly part to be pressurized and compression speed controller.

Touch key pad display.

All human body is supplied by peripheral circulation of blood and lymph. Our systems help patients to recover from their disease by improving blood and lymph circulation and metabolism. Our systems consist of air chambers and an air pump offering air compression massage to legs, arms or waist.

Medical indications

- DVT (Deep Vein Thrombosis)
- Lymph edema
- Post-paralytic
- Limb convulsion
- Venosclerosis
- Chronic venous insufficiency (CVI)
- Post-mastectomy
- Post-trauma
- Spinal cord injury
- Varicose veins
- Sciatica
- Sprains
- Diabetes
- Arterial insufficiency

Home care indications

- Lymph edema
- Diseases caused by interrupted blood circulation
- Cellulitis fat dissolution
- Pregnant woman's limb treatment
- Numb feet or hands
- Cold feet or hands
- Settlement of stress
- Exercise effects
- Insomnia
- Healthier body and elastic skin
- Old age, weakness and fatigue

Manual: GB, IT, ES

Technical Specifications

- Pressure range: 10~200 mmHg
- Automatic pressure detection
- Operating voltage: 220-240 V, 50/60 Hz, 50 VA
- Use time: 1~90 minutes
- Body size: 415 x 310 x 160 mm
- Body weight: 6.5 kg