

# 3-SECTION PILLOW - 100% cotton



Product Code	28551
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

EAN/UPC	8053470771972
PARAF	975385941

## Description

### 3-SECTION CUSHION

Size: 45 x 45 cm

Products to prevent low decubitus risk, made of silicone fibre, coated with 100% preshrunk cotton cover. Main property of silicone fibre is to maintain original shape also after strong pressure. This allows to maintain comfort and proper weight distribution in the time. Breathable cotton cover provides maximum cool comfort. Washable at 30-40°C.

All cushions and pillow are supplied with multilingual instructions: GB, FR, IT, ES, PT, DE

Made in Italy.