

# PILLOW WITH CONCAVE CENTRE - 100% cotton



Product Code	28552
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

PARAF	975385954
-------	-----------

## Description

### CUSHION WITH CONCAVE CENTRE

Size: 45 x 45 cm

Products to prevent low decubitus risk, made of silicone fibre, coated with 100% preshrunk cotton cover. Main property of silicone fibre is to maintain original shape also after strong pressure. This allows to maintain comfort and proper weight distribution in the time. Breathable cotton cover provides maximum cool comfort. Washable at 30-40°C.

All cushions and pillow are supplied with multilingual instructions: GB, FR, IT, ES, PT, DE

Made in Italy.